

Aqua Stim™

Water Caloric Irrigator

Compact and quiet, the Micromedical designed water caloric irrigator has an easy to fill and clean external water tank providing water for your next irrigation. VisualEyes™ VNG software automatically selects irrigation temperature for your chosen test. The irrigator handle pushbutton simultaneously starts nystagmus recording and water flow.



- ▶ Stand-alone or PC USB Interface
- ▶ Micromedical's Spectrum software control
- ▶ Time to temperature: less than 3 minutes
- ▶ Handle features LED's to illuminate ear; pushbutton switch; disposable tips
- ▶ External water tank



Features

Water Flow Rate: 250cc/30 seconds
Water Temperatures: 30°C and 44°C
Water fill and heat cycle: Each irrigation
Water hose length: 3 m (9.8 ft)

Temperature Stability: ± 1°C
Irrigation time: 30 seconds
Flow Accuracy: ± 15 ml/irrigation
Accuracy at Tip: ± 1°C

Specifications

Case Dimensions: (WxDxH) 35cm x 30cm x 21cm
 13.8in x 12in x 8.3in
Voltage: 110-130 VAC @ 50-60 Hz
 220-240 VAC @ 50-60 Hz
Power: 600 VA

External water tank: 3.5 liters / 0.9 gallons
Weight of irrigator: 5.8 kg / 12.8 lbs
Chassis Leakage: < 300 µA (single fault)