

Handheld Screening Tympanometer

by Maico Diagnostics

easyTymp[®]

Introducing the easyTymp[®]

Handheld tympanometer with a unique, space saving design. Fast, accurate, middle ear measurements.

easyTymp[®] without printer and removable ear tip box pictured



*The ideal handheld device for:
Audiology, Family Medicine,
Pediatrics, and School settings.*

- NOAH Compatible
- Ergonomic Design
- Full color; high resolution graphic result display
- Probe Tone; 226 Hz (standard)
– Optional; 1000 Hz
- Ipsi Reflex Frequencies; 500 HZ, 1KHz, 2KHz, 4KHz
- Rechargeable docking station; supplied with extra battery

www.maico-diagnostics.com



easyTymp®

The easyTymp® is designed to make detection and documentation of middle ear pathologies fast and accurate.



easyTymp® with printer pictured



easyTymp® carrying case pictured

Equipment Configurations:

- easyTymp® with printer #721B
- easyTymp® without printer #721C

Optional Equipment:

- Wall mount #6699
- Carrying case #1035-3024
- 1000 Hz Probe tone #721-14

Accessories:

- Paper roll #5529
- Eartip kit #721-5
- Rechargeable battery #721-4

Impedance Measuring System

Probe Tone	Frequency:	226 Hz, 1000 Hz (optional)
Air Pressure	Control:	Automatic
	Range:	-300 to +200 daPa
	Pressure limitation:	-750 daPa and +550 daPa
Compliance	Range:	0.1 to 8.0 ml at 226 Hz probe tone (Ear volume: 0.1 to 8.0 ml) and 0.1 to 15 mmho at 1000 Hz probe tone
Test Types	Tympanometry:	Automatic, where the start and stop pressure can be user-programmed in the set-up function

Reflex Functions

Signal Sources	Tone-Ipsi, Reflex:	500, 1000, 2000, 4000 Hz
	Noise-Ipsi, Reflex:	Wide band
Test Types	Automated Reflex:	Automatic reflexes: - Single screening intensity or Reflex auto search

This brochure contains only a small segment of the comprehensive product portfolio of MAICO. To find out more about other solutions, please contact us at: Telephone: 888-941-4201



MAICO Diagnostics
7625 Golden Triangle Drive
Eden Prairie, MN 55344
Tel: 888-941-4201
Fax: 952-903-4100
Web: www.maico-diagnostics.com

ONE YEAR LIMITED WARRANTY: This warranty is extended to the original purchaser of the instrument, by Maico, through the Distributor from whom it was purchased. The warranty covers defects in material and workmanship for a period of one year from date of invoice of the instrument to the original purchaser. Accessories which are purchased from Maico at the same time as the instrument are warranted for one year from the date of invoice. For additional information contact your Maico Distributor.