

MA 1

Air Conduction Audiometer



The super-lightweight MA1 weighs in at only 1 pound. Small enough to hold in the palm of your hand! Tucked into its carrying case (included), the MA1 is our smallest, lightest portable audiometer.

- **Air conduction screener**
- **Easy one hand push button operation**
- **Traditional binaural TDH-39 headset included**
- **Fits in the palm of your hand**
- **Exceptionally lightweight, only 1 pound**
- **Carrying case included**

EXCELLENCE IN HEARING SCREENING SINCE 1937



MA 1

Accurate screening in the palm of your hand

SLEEK, SMALL DESIGN

With its super-lightweight design, the MA 1 is small enough to hold in the palm of your hand. Tucked into its carrying case (included) the MA1 is our smallest, lightest portable audiometer.

ACCURATE, DEPENDABLE ENGINEERING

Because it's a portable, the MA 1 is engineered to be tough and reliable. Take the MA 1 anywhere with confidence. It's user-friendly push button design makes adjusting frequency and intensity easy.

THE NECESSITIES OF HEARING SCREENING

Although the MA 1 is small, it offers a full range of testing screening. Test frequencies 500, 1000, 2000 and 4000Hz. Signal intensities range from 15 to 50 dB HL.

GAIN NEW TESTING FLEXIBILITY

To provide you with more versatile testing opportunities, the MA 1 is powered by 2 AA alkaline batteries, only.

* 2.1 lb with carrying case



TECHNICAL DATA

Frequency Range_____500, 1000, 2000, 4000Hz

Intensity Range_____15 to 50 dB HL

Calibrated to ANSI S3.6

CPT CODES

92551_____Screening, pure tone

92552_____Threshold, air

ONE YEAR LIMITED WARRANTY: This warranty is extended to the original purchaser of the instrument, by Maico, through the Distributor from whom it was purchased. The warranty covers defects in material and workmanship for a period of one year from date of delivery of the instrument to the original purchaser. Accessories which are purchased from Maico at the same time as the instrument are warranted for one year from the date of purchase. For additional information contact your Maico Distributor.

EXCELLENCE IN HEARING SCREENING SINCE 1937

