

Race Car[®] Tympanometer

Successfully Screen Young Children in Seconds



*The Maico Race Car[®]
Tympanometer*

The only tympanometer specifically designed to capture and keep your young patient's attention.



- **Visual "Race Car" distraction and multimedia entertainment holds the child's attention during the test**
- **Choose the Tympanometer or the Tympanometer/Audiometer combination**
- **Perform a tympanogram with ipsilateral reflex**
- **Quick, quiet internal printer**



EXCELLENCE IN HEARING SCREENING SINCE 1937

Race Car®

Quick, Accurate, Professional & Fun

THE RACE CAR TYMPANOMETER OR THE RACE CAR TYMPANOMETER/AUDIOMETER

Hearing screening can be difficult. A child who is uncooperative can be a time consuming challenge. The Race Car Tympanometer provides a visual distraction that captures the child's attention so you can create and maintain the proper seal for an accurate tympanogram. The entire test takes place while the patient watches the car accelerate across the screen to the finish line.

An accurate tympanogram with ipsilateral reflex is ready to review or print in seconds. The Race Car Tympanometer is also available as a Tympanometer/Audiometer combo.



TECHNICAL DATA

TYMPANOMETRY MODE

Probe Frequency _____ 226 Hz, 1000 Hz (optional) \pm 3%

REFLEX MODE

Reflex Tones _____ 500, 1000, 2000, 4000 Hz

AUDIOMETRY MODE

Frequencies _____ 125 to 8000 Hz

Level _____ -10 to 90 dBHL

Specifications subject to change without notice



Visual Race Car and Multimedia Software Program

EXCELLENCE IN HEARING SCREENING SINCE 1937

