



Monitor Instruments, Inc. is a manufacturer of industrial audiometers and audiometric accessories, offering quality instruments and professional calibration services for audiometers and acoustic safety equipment since the early 1970's.

(800) 853-6785 · (919) 732-5400 · FAX (919) 732-6153 · info@monitorinstrumentsinc.com

The MI-300

The MI-300 Calibration Monitor performs a quick and efficient daily biological calibration check of your audiometer, as required by OSHA. Its lightweight and durable construction will insure years of trouble-free operation.



Features

- LED indicators for right response, left response, and power on.
- Calibration at 70 dB hearing level at 1 kHz.
- Both normally open and normally closed response relay configuration.
- Automatic power off when not in use.
- Designed for wall mount or table top use.
- Works with manual or automatic microprocessor audiometers.

General Description

- Power: 9 volt alkaline battery
- Size: 6.0 x 6.0 x 3.3 inches
- Weight: 1.0 lb.
- Construction: Aluminum housing with delrin coupler.

Standard Accessories

- 9 volt battery
- Response cable (specify make and model of audiometer when ordering)
- Rubber feet
- Operator manual