

The REINFORCER **WR-400 Wired VRA** system combines motion and attention getting light patterns to successfully perform pediatric hearing testing. Each smoked Plexiglass box features a friendly cartoon character that lights up and moves back and forth rapidly and quietly when activated. Additionally a horseshoe shaped light bar with 20 bright LEDs flash sequentially or in a marquee pattern to vary the combination of lights and motion reinforcement. The lightweight Plexiglass boxes are triangular shaped to mount easily in the corners of the sound room with Velcro or double back tape. This user installable system operates at low voltages for safety and uses modular (telephone style) plug in connectors to easily install the system into a sound room. Optional centering box and foot switch are available.



KEY FEATURES

- EASY TO INSTALL
- TRIANGULAR SHAPED BOXES FIT COMPLETELY INTO CORNERS
- WAVE/MARQUEE LIGHT PATTERNS
- FULL MOTION ACTION
- SEPERATE MOTOR (ACTION) CONTROL -NEW

SYSTEM INCLUDES

- Two Triangular Plexiglass Boxes
- Wireless Hand-Held Multifunction Remote Control
- Plug-In Wall Transformer
- 20 ft Long Interface Cables

OPTIONAL ACCESSORIES

- Centering Box
- Installation Kit
- Foot Switch



WR-400 Wired Control Console

